

Vivago FIDO™ Personal alarm device

Quick reference guide

If you feel that you are in need of help, you should:

- Immediately press the alarm button on the device
- Press until the red alarm LED turns on

Important!

- The device notifies the user about battery charge getting low by blinking the red LED once in every 16 seconds.
- The device also notifies the system when the battery must be charged or when it is full.
- After the battery low notification, the battery must be charged within three days.
- Alarms can be sent at the coverage area of the Vivago institution base unit. Thick stonewalls and metal constructions may prevent data transfer between the device and the base unit.
- The device is waterproof, so you can keep it with you in the shower or in the bath. The device can't be used in hot conditions, such as Sauna.
- Familiarize yourself with the Vivago FIDO™ User Manual before using the device.



Vivago

Vivago Ltd

Itsehallintokuja 6, FI-02600 Espoo, Finland
Tel. +358 10 2190 610 | email: info@vivago.com | www.vivago.com